



AEL IR 4400 – The World’s Best Equine Laser

If you’ve ever wondered why some horses consistently win and perform better than others? The answer is simple. They are being treated with Class 3B laser therapy!

The AEL IR 4400 is an Australian designed, engineered and manufactured laser that is used to treat all common equine injuries.

Australian Equine Lasers developed the world’s best equine laser, the AEL IR 4400, with a 904nm Cluster head and a massive 450mW of output power. The 4400 uses advanced laser diode technology to effectively deliver laser dosages that heal injuries.

How does the 4400 work?

The AEL IR 4400 uses the best wavelength, 904nm to treat a wide variety of ailments. Using the 4400, you can achieve drug like results without the harmful side effects associated with taking prescription medications.

The AEL IR 4400 delivers the required energy at the injured muscle group, encouraging blood flow, reducing pain, inflammation and healing damaged cells. Recovery times are accelerated.

The 4400 is successfully used to treat:

- Arthritis
- Back pain
- Injured tendons, muscles, ligaments
- Fetlocks
- Pastern injuries
- Wounds
- Joints
- Scars
- As well as many more!

Best of all, the 4400 requires no formal training!

Need further proof? Read the 5 star reviews and glowing customer testimonials on Facebook from people who are successfully using the 4400 to treat their injured horse.

https://www.facebook.com/Australian-Equine-Lasers-122628175805396/reviews/?ref=page_internal

If you’re interested in learning more or to place your order contact Tracey at Australian Equine Lasers on **0413 936 976**.

Winning Is Everything